epicite^{hydro}

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The Ideal Burn Dressing?

92% of respondents:

- there is no ideal burn dressing available

Most important properties:

- Pain-free (100%)
- Non-adhesion (87.5%)
- Absorbency (87.5%)
- Antimicrobial activity (83.5%)

BURNS 38 (2012) 960-966

Available online at www.sciencedirect.com

SciVerse ScienceDirect

journal homepage: www.elsevier.com/locate/burns

The properties of an “ideal” burn wound dressing – What do we need in daily clinical practice? Results of a worldwide online survey among burn care specialists

Harald F. Selig^{a,b}, David B. Lumenta^{c,d}, Michael Giretzlehner^e,
Marc G. Jeschke^f, Dominic Upton^g, Lars P. Kamolz^{a,c,*}

^a Section of Plastic, Aesthetic and Reconstructive Surgery, Department of Surgery, Landeskrankenhaus Wiener Neustadt, Austria
^b Hospital for Hand Surgery Bad Neustadt, Bad Neustadt/Saale, Germany
^c Division of Plastic, Aesthetic and Reconstructive Surgery, Department of Surgery, Medical University of Graz, Austria
^d Division of Plastic and Reconstructive Surgery, Department of Surgery, Medical University of Vienna, Austria
^e Research Unit for Medical-Informatics, RISC Software GmbH, Johannes Kepler University Linz, Linz, Austria
^f Ross Tilley Burn Centre, Sunnybrook Health Sciences Centre, Sunnybrook Research Institute, Division of Plastic Surgery, Department of Surgery, University of Toronto, Toronto, Canada
^g Institute of Health and Society, University of Worcester, England, United Kingdom

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epicite^{hydro} zur Wundbett-Konditionierung

1 – 2a°
superficial to
superficial
partial thickness

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




epicite^{hydro} 10 x 10 cm


Produkt A

SEM = Scanning Electron Microscopy at magnitude of 3,00 K X

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

Figure 1. Uptake of PHMB (Prontosan) into epicite^{hydro}

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Volume uptake

epicite ^{hydro} size [cm]	volume of liquid in epicite ^{hydro} [ml]	
	from package	after 2h incubation
10 x 10	19	26
15 x 20	56	78
20 x 20	74	104

Uptake of antiseptic compound

Antisepticum	concentration of compound [mg/L]				Minimal inhibition conc. vs. MRSA [mg/L]
	original solution	in epicite ^{hydro} , @ 30 min	in epicite ^{hydro} , @ 60 min	in epicite ^{hydro} , @ 120 min	
Lavanid 2 [PHMB]	400	110	190	230	0,5-1,0 ^b
Prontosan [PHMB]	1.000	260	530	760	
Octenisept [Octinidine]	1.000	360	610	800	0,375 - 3,0 ^c
Polysept [Povidone iodine]	100.000	10.000	20.000	50.000	0,25 - 1,0 ^d
Mafenide [Mafenide]	50.000	20.000	33.000	40.000	1 ^e

Table 1. Uptake of different antiseptics into epicite^{hydro}, size, mean liquid volumes from package (0,19 ml/cm²) and after 2 h incubation (0,26 ml/cm²). Concentrations of drug substances in the original solution and after uptake from five volumes (200 ml) for 30 min, 60 min and 120 min compared to the minimal inhibition concentrations against MRSA.

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QRSKIN GmbH